



GE Medical Systems

GE HISPEED ADVANTAGE CT/I CT SCANNER

GENERAL NOTES

- ALL SPECIFICATIONS AND SPACE REQUIREMENTS ARE SHOWN GENERICALLY. CONTACT GAMMA TECH, INC. FOR SPECIFIC REQUIREMENTS PERTAINING TO OPTIONAL EQUIPMENT COMPONENTS AND ADDITIONAL ROOM PREPARATION SPECIFICATIONS.
- RECOMMENDED FINISHED CEILING HEIGHT IS 9'-0". MINIMUM FINISHED CEILING HEIGHT 8'-6".
- AMBIENT OPERATION TEMPERATURE: 55° TO 75° F. HUMIDITY: 20% TO 80% (NON-CONDENSING). ALLOWABLE TEMPERATURE CHANGE: 15°F PER HOUR. ALLOWABLE HUMIDITY CHANGE: 10% PER HOUR. MAX - ALTITUDE: 8000FT ABOVE SEA LEVEL.
- MAIN INCOMING POWER SUPPLY: 480VAC, THREE PHASE, 112.5KVA, DELTA OR WYE, 90AMP, 3 POLE, MAIN CB 5% VOLTAGE REGULATION, 95% POWER FACTOR

© PLANNING DIVISION


FILE NAME:
GE HISPEED ADVANTAGE
CT/I

REVISION DATE:
01.22.2008


EQUIPMENT SCHEDULE

EQUIPMENT IDENTIFICATION NUMBER

ID	BTU / HOUR	DESCRIPTION	WEIGHT (LBS.)			
			DIMENSIONS (IN.)			
			W	x	D	x
1	24403	GENERAL ELECTRIC MEDICAL SYSTEMS HISPEED ADVANTAGE CT/I CT SCANNER GANTRY	88.00	42.75	78.55	3700
2	700	GENERAL ELECTRIC MEDICAL SYSTEMS HISPEED ADVANTAGE CT/I CT SCANNER PATIENT TABLE W/ EXTENDED TABLE	24.25	91.91	43.78	1185
3	4500	GENERAL ELECTRIC MEDICAL SYSTEMS HISPEED ADVANTAGE CT/I CT SCANNER OPERATOR CONSOLE	48.00	40.00	48.00	662
4	3400	GENERAL ELECTRIC MEDICAL SYSTEMS HISPEED ADVANTAGE CT/I SCANNER POWER DISTRIBUTION UNIT	33.00	18.50	78.00	800
5	N/A	MEDRAD STELLANT CT INJECTOR SYSTEM HEAD ON CEILING MOUNT COUNTERPOISE				
6	N/A	MEDRAD STELLANT CT INJECTOR SYSTEM ELECTRONICS BASE UNIT (ON COUNTERTOP)	11.00	8.80	11.50	14
7	N/A	MEDRAD STELLANT CT INJECTOR SYSTEM DISPLAY CONTROL CONSOLE (ON COUNTERTOP)	12.50	8.90	13.50	8
8	406	GE MEDICAL SYSTEMS ADVANTAGE WINDOWS STATION				105



TYPICAL ROOM AND EQUIPMENT LAYOUT
(DASHED LINES DEPICT EQUIPMENT TRAVEL AND SERVICE CLEARANCES)
RECOMMENDED FINISHED CEILING HEIGHT = 9'-0"



PLANNING BY:



GAMMA TECH, INC

129 Wild Berry Lane
Hampstead, NC 28443
Ph (910) 270-9787 • Fax (910) 270-9767
www.gammatech.net