





GE Medical Systems

HISPEED NXi CT SCANNER



TYPICAL ROOM AND EQUIPMENT LAYOUT
(DASHED LINES DEPICT EQUIPMENT TRAVEL AND SERVICE CLEARANCES)
RECOMMENDED FINISHED CEILING HEIGHT IS 9'-0"



GENERAL NOTES

- ALL SPECIFICATIONS AND SPACE REQUIREMENTS ARE SHOWN GENERICALLY. CONTACT GAMMA TECH, INC. FOR SPECIFIC REQUIREMENTS PERTAINING TO OPTIONAL EQUIPMENT COMPONENTS AND ADDITIONAL ROOM PREPARATION SPECIFICATIONS.
- RECOMMENDED FINISHED CEILING HEIGHT IS 9'-0".
MINIMUM FINISHED CEILING HEIGHT IS 7'-6".
- AMBIENT OPERATION TEMPERATURE: 68° TO 82° F.
HUMIDITY: 30% TO 70% (NON-CONDENSING).
ALLOWABLE TEMPERATURE CHANGE: 8°F PER HOUR
ALLOWABLE HUMIDITY CHANGE: 5% PER HOUR
MAX - ALTITUDE: 8,200FT ABOVE SEA LEVEL
- MAIN INCOMING POWER SUPPLY:
480VAC, THREE PHASE, 112.5KVA, DELTA OR WYE, 90AMP, 3 POLE, MAIN CB
5% VOLTAGE REGULATION, 90% POWER FACTOR

© PLANNING DIVISION

FILE NAME:
GE HISPEED NXi

REVISION DATE:
07.02.2008

EQUIPMENT SCHEDULE

EQUIPMENT IDENTIFICATION NUMBER

EQUIPMENT IDENTIFICATION NUMBER	BTU / HOUR	DESCRIPTION	WEIGHT (LBS.)			
			DIMENSIONS (IN.)			
			W	x	D	x
1	16260	GENERAL ELECTRIC MEDICAL SYSTEMS HISPEED NXi GANTRY	71.65	35.32	72.84	2580
2	318	GENERAL ELECTRIC MEDICAL SYSTEMS HISPEED NXi PATIENT TABLE WITHOUT EXTENDED TABLE TOP	25.60	88.32	37.40	651
3	2911	GENERAL ELECTRIC MEDICAL SYSTEMS HISPEED NXi OPERATOR'S CONSOLE / COMPUTER	43.31	33.90	51.90	449
4	615	GENERAL ELECTRIC MEDICAL SYSTEMS HISPEED NXi POWER DISTRIBUTION UNIT	21.66	27.56	32.29	346

PLANNING BY:



GAMMA TECH, INC

129 Wild Berry Lane
Hampstead, NC 28443
Ph (910) 270-9787 • Fax (910) 270-9767
www.gammatech.net