





GE Medical Systems

GE LIGHTSPEED RT 16 CT SCANNER WITH GT1700 TABLE



TYPICAL ROOM AND EQUIPMENT LAYOUT
(DASHED LINES DEPICT EQUIPMENT TRAVEL AND SERVICE CLEARANCES)
RECOMMENDED FINISHED CEILING HEIGHT IS 9'-0"



GENERAL NOTES

- ALL SPECIFICATIONS AND SPACE REQUIREMENTS ARE SHOWN GENERICALLY. CONTACT GAMMA TECH, INC. FOR SPECIFIC REQUIREMENTS PERTAINING TO OPTIONAL EQUIPMENT COMPONENTS AND ADDITIONAL ROOM PREPARATION SPECIFICATIONS.
- RECOMMENDED FINISHED CEILING HEIGHT IS 9'-0".
MINIMUM FINISHED CEILING HEIGHT IS 8'-0".
- AMBIENT OPERATION TEMPERATURE: 64° TO 79° F.
HUMIDITY: 30% TO 60% (NON-CONDENSING).
ALLOWABLE TEMPERATURE CHANGE: 5°F PER HOUR
ALLOWABLE HUMIDITY CHANGE: 5% PER HOUR
MAX - ALTITUDE: 10,000FT ABOVE SEA LEVEL
- MAIN INCOMING POWER SUPPLY:
480VAC, THREE PHASE, 225KVA, WYE CONNECTED, 125AMP, 3 POLE, MAIN CB
5% VOLTAGE REGULATION, 85% POWER FACTOR

© PLANNING DIVISION

FILE NAME:
GE LIGHTSPEED RT 16

REVISION DATE:
07.02.2008

EQUIPMENT SCHEDULE

EQUIPMENT IDENTIFICATION NUMBER

ID	BTU / HOUR	DESCRIPTION	WEIGHT (LBS.)					
			DIMENSIONS (IN.)					
			W	x	D	x	H	
1	43000	GENERAL ELECTRIC MEDICAL SYSTEMS LIGHTSPEED RT 16 CT SCANNER GANTRY	96.00		39.60		77.40	3891
2	700	GENERAL ELECTRIC MEDICAL SYSTEMS LIGHTSPEED RT 16 CT SCANNER GT 1700 PATIENT TABLE WITH EXTENDED TABLE TOP	25.80		92.30		43.10	881
3	8631	GENERAL ELECTRIC MEDICAL SYSTEMS LIGHTSPEED RT 16 CT SCANNER OPERATOR'S CONSOLE / COMPUTER	48.73		48.35		34.76	493
4	3400	GENERAL ELECTRIC MEDICAL SYSTEMS LIGHTSPEED RT 16 CT SCANNER POWER DISTRIBUTION UNIT	27.60		21.70		41.70	701
5	1109	GENERAL ELECTRIC MEDICAL SYSTEMS ADVANTAGESIM MD - VIRTUAL SIMULATION WORKSTATION WITH TWO LCD MONITORS						81
6	N/A	MEDRAD STELLANT CT INJECTOR SYSTEM INJECTOR HEAD ON STAND	6.10		14.50		65.30	112
7	N/A	MEDRAD STELLANT CT INJECTOR SYSTEM ELECTRONICS BASE UNIT (COUNTER MOUNTED)	11.00		8.80		11.50	14
8	N/A	MEDRAD STELLANT CT INJECTOR SYSTEM DISPLAY CONTROL CONSOLE (COUNTER MOUNTED)	12.50		8.90		13.50	8
9		LAP OF AMERICA MODEL CT1-1 CEILING MOUNTED MOVABLE PATIENT POSITION LASER	43.37		7.06		7.00	76
10		LAP OF AMERICA MODEL CT1-1 WALL MOUNTED FIXED PATIENT POSITION LASER	4.33		3.98		8.11	5
11		LAP OF AMERICA MODEL CT1-1 KEYPAD	5.00		1.50		5.50	1

PLANNING BY:



GAMMA TECH, INC

129 Wild Berry Lane
Hampstead, NC 28443
Ph (910) 270-9787 • Fax (910) 270-9767
www.gammatech.net