





# GE Medical Systems

## GE PROSPEED PLUS CT SCANNER



**TYPICAL ROOM AND EQUIPMENT LAYOUT**  
(DASHED LINES DEPICT EQUIPMENT TRAVEL AND SERVICE CLEARANCES)  
RECOMMENDED FINISHED CEILING HEIGHT = 8'-0"



### GENERAL NOTES

- ALL SPECIFICATIONS AND SPACE REQUIREMENTS ARE SHOWN GENERICALLY. CONTACT GAMMA TECH, INC. FOR SPECIFIC REQUIREMENTS PERTAINING TO OPTIONAL EQUIPMENT COMPONENTS AND ADDITIONAL ROOM PREPARATION SPECIFICATIONS.
- RECOMMENDED FINISHED CEILING HEIGHT IS 8'-0". MINIMUM FINISHED CEILING HEIGHT 7'-6".
- AMBIENT OPERATION TEMPERATURE: 60° TO 75° F. HUMIDITY: 30% TO 60% (NON-CONDENSING). ALLOWABLE TEMPERATURE CHANGE: 5°F PER HOUR. ALLOWABLE HUMIDITY CHANGE: 10% PER HOUR. MAX - ALTITUDE: 10,000FT ABOVE SEA LEVEL
- MAIN INCOMING POWER SUPPLY:  
480VAC, THREE PHASE, 112.5KVA, DELTA OR WYE, 90AMP, 3 POLE, MAIN CB  
5% VOLTAGE REGULATION, 90% POWER FACTOR

© PLANNING DIVISION

FILE NAME:  
GE PROSPEED PLUS

REVISION DATE:  
01.22.2008

## EQUIPMENT SCHEDULE

### EQUIPMENT IDENTIFICATION NUMBER

ID	DESCRIPTION	BTU / HOUR		WEIGHT (LBS.)	
		DESCRIPTION		DIMENSIONS (IN.)	
		W	x	D	x
1	GENERAL ELECTRIC MEDICAL SYSTEMS PROSPEED SERIES CT SCANNER GANTRY	84.65	37.40	74.41	2536
2	GENERAL ELECTRIC MEDICAL SYSTEMS PROSPEED SERIES CT SCANNER PATIENT TABLE	103.54	22.44	40.00	948
3	GENERAL ELECTRIC MEDICAL SYSTEMS PROSPEED SERIES CT SCANNER OPERATOR CONSOLE (Zj-a)	48.82	32.28	50.00	507
4	GENERAL ELECTRIC MEDICAL SYSTEMS PROSPEED SERIES CT SCANNER POWER DISTRIBUTION UNIT	21.61	24.80	32.32	673
5	MEDRAD STELLANT CT INJECTOR SYSTEM HEAD ON CEILING MOUNT COUNTERPOISE				
6	MEDRAD STELLANT CT INJECTOR SYSTEM ELECTRONICS BASE UNIT (ON DESK)	11.00	8.80	11.50	14
7	MEDRAD STELLANT CT INJECTOR SYSTEM DISPLAY CONTROL CONSOLE (ON DESK)	12.50	8.90	13.50	8

PLANNING BY:



**GAMMA TECH, INC**

129 Wild Berry Lane  
Hampstead, NC 28443  
Ph (910) 270-9787 • Fax (910) 270-9767  
www.gammatech.net