




GE Medical Systems

GE PROSPEED S/SX CT SCANNER



TYPICAL ROOM AND EQUIPMENT LAYOUT
(DASHED LINES DEPICT EQUIPMENT TRAVEL AND SERVICE CLEARANCES)
RECOMMENDED FINISHED CEILING HEIGHT = 8'-0"

GENERAL NOTES

- ALL SPECIFICATIONS AND SPACE REQUIREMENTS ARE SHOWN GENERICALLY. CONTACT GAMMA TECH, INC. FOR SPECIFIC REQUIREMENTS PERTAINING TO OPTIONAL EQUIPMENT COMPONENTS AND ADDITIONAL ROOM PREPARATION SPECIFICATIONS.
- RECOMMENDED FINISHED CEILING HEIGHT IS 8'-0". MINIMUM FINISHED CEILING HEIGHT 7'-6".
- AMBIENT OPERATION TEMPERATURE: 60° TO 75° F. HUMIDITY: 30% TO 60% (NON-CONDENSING). ALLOWABLE TEMPERATURE CHANGE: 5°F PER HOUR. ALLOWABLE HUMIDITY CHANGE: 10% PER HOUR.
MAX - ALTITUDE: 10,000FT ABOVE SEA LEVEL
- MAIN INCOMING POWER SUPPLY:
480VAC, THREE PHASE, 112.5KVA, DELTA OR WYE, 90AMP, 3 POLE, MAIN CB
5% VOLTAGE REGULATION, 95% POWER FACTOR

© PLANNING DIVISION

FILE NAME:
GE PROSPEED S-SX

REVISION DATE:
01.24.2008

EQUIPMENT SCHEDULE

EQUIPMENT IDENTIFICATION NUMBER

EQUIPMENT IDENTIFICATION NUMBER	BTU / HOUR	DESCRIPTION	WEIGHT (LBS.)			
			DIMENSIONS (IN.)			
			W	x	D	x
1	11150	GENERAL ELECTRIC MEDICAL SYSTEMS PROSPEED SX CT SCANNER GANTRY	84.65	37.40	74.41	2536
2	400	GENERAL ELECTRIC MEDICAL SYSTEMS PROSPEED SX CT SCANNER PATIENT TABLE W/ EXTENDED TABLE	103.54	22.44	40.00	948
3	1985	GENERAL ELECTRIC MEDICAL SYSTEMS PROSPEED SX CT SCANNER OPERATOR CONSOLE (Qj-b)	31.50	29.72	47.36	250
4	3584	GEMS PROSPEED SX CT SCANNER POWER DISTRIBUTION UNIT (MUST BE LOCATED A MIN. OF 6'-0" FROM GANTRY AND/OR CONSOLE)	21.61	24.80	32.32	673
5	N/A	MEDRAD STELLANT CT INJECTOR SYSTEM HEAD ON CEILING MOUNT COUNTERPOISE				
6	N/A	MEDRAD STELLANT CT INJECTOR SYSTEM ELECTRONICS BASE UNIT (ON DESK)	11.00	8.80	11.50	14
7	N/A	MEDRAD STELLANT CT INJECTOR SYSTEM DISPLAY CONTROL CONSOLE (ON DESK)	12.50	8.90	13.50	8

PLANNING BY:



GAMMA TECH, INC

129 Wild Berry Lane
Hampstead, NC 28443
Ph (910) 270-9787 • Fax (910) 270-9767
www.gammatech.net

