



PLANNING BY:



CAMMA TECH, INC

Hampstead, NC 28443 Ph (910) 270-9787 • Fax (910) 270-9767 www.gammatech.net



GENERAL NOTES

ALL SPECIFICATIONS AND SPACE REQUIREMENTS
ARE SHOWN GENERICALLY. CONTACT GAMMA
TECHNOLOGIES FOR SPECIFIC REQUIREMENTS
PERTAINING TO OPTIONAL EQUIPMENT COMPONENTS
AND ADDITIONAL ROOM PREPARATION SPECIFICATIONS.

© PLANNING DIVISION FILE NAME: GE CPI RAD REVISION DATE: 01.31.2008

EQUIPMENT SCHEDULE

EQUIPMENT IDENTIFICATION NUMBER

			BTU / HOUR		WEIGHT (LBS.)		
NOTES			DESCRIPTION	DIN w	MENSIOI × D	NS (IN.)	
OFOI	1	150	CPI INDICO 100 SERIES HIGH FREQUENCY GENERATOR CONTROL CONSOLE (COUNTER MOUNTED)	16.60	12.30 NOT INCLUDE	3.70	11
OFOI	2	4000	CPI INDICO 100 SERIES HIGH FREQUENCY GENERATOR CONTROL CABINET	18.00	16.00	48.70	235
OFOI	3	1200	GENERAL ELECTRIC MEDICAL SYSTEMS RFX 90/90 R&F TABLE	90.00	68.00	60.00	3000
OFOI	4	1200	GENERAL ELECTRIC MEDICAL SYSTEMS RFX TABLE UNIVERSAL SYSTEMS CABINET	33.00	15.00	79.00	950
OFOI	5	1200	GENERAL ELECTRIC MEDICAL SYSTEMS XT RADIOGRAPHIC SUSPENSION WITH INBOARD MOUNTINGS	228.00 N	120.50 MODIFIED TO	94.00 FIT ROOM BY SUI	941 PPLIER
OFOI	6	200	GENERAL ELECTRIC MEDICAL SYSTEMS SENTRY III COLLIMATOR ELECTRONICS CABINET (WALL MOUNTED)	17.00	8.00	26.00	50
OFOI	7	200	PAUSCH BS-2000 V VERTICAL BUCKY STAND	31.61	16.73	87.79	300
OFOI	8		H.S. LIONS MODEL610 MONITOR CART WITH KEYBOARD SHELF CONTAINING THE FOLLOWING:	34.00	34.00	57.00	175
OFOI		300	INFIMED PLATINUM ONE 21" IMAGE MONITOR	16.00	19.00	17.00	
OFOI			INFIMED PLATINUM ONE PC KEYBOARD	19.00	8.00	3.00	
OFOI	9	800	INFIMED PLATINUM ONE (BENEATH COUNTER) PC CABINET	7.75	18.13	16.50	29
OFOI	10	140	INFIMED PLATINUM ONE (BENEATH COUNTER) ISOLATION TRANSFORMER	8.50	12.25	7.25	32
OFOI	11	140	INFIMED PLATINUM ONE (BENEATH COUNTER) GENERATOR INTERFACE MODULE	2.56	12.50	12.50	6
OFOI	12	140	INFIMED PLATINUM ONE (BENEATH COUNTER) UNINTERRUPTED POWER SUPPLY	7.00	17.10	8.80	57

NOTE EXPLANATION: CFCI = CONTRACTOR FURNISHED - CONTRACTOR INSTALLED OFCI = OWNER FURNISHED - CONTRACTOR INSTALLED OFOI = OWNER FURNISHED - OWNER INSTALLED