



TYPICAL ROOM LAYOUT



GENERAL NOTES

- ALL SPECIFICATIONS AND SPACE REQUIREMENTS ARE SHOWN GENERICALLY. CONTACT GAMMA TECH, INC. FOR SPECIFIC REQUIREMENTS PERTAINING TO OPTIONAL EQUIPMENT COMPONENTS AND ADDITIONAL ROOM PREPARATION SPECIFICATIONS.
- MINIMUM RECOMMENDED CEILING HEIGHT IS 8'-6".
- SUGGESTED AMBIENT ROOM REQUIREMENTS ARE 66° TO 75° F. AT 30% TO 60% R.M. (NON-CONDENSING).
- MAIN INCOMING POWER SUPPLY: VARIES WITH GENERATOR MODEL.

GENDEX HIGH FREQUENCY RADIOGRAPHIC SYSTEM SUITE.

EQUIPMENT SCHEDULE

		EQUIPMENT IDENTIFICATION NUMBER				
		BTU / HOUR	WEIGHT (LBS.)			
		DESCRIPTION	DIMENSIONS (IN.)			
			W	D	H	
1	320	GENDEX ATC-SERIES HF GENERATOR OPERATOR CONTROL PANEL	16.00	12.75	36.50	80
2	4000	GENDEX ATC-SERIES HF GENERATOR RACK CABINET WITH H.T. TRANSFORMER	26.50	19.50	63.25	250
3	393	GENDEX EV-200 ELEVATING TABLE	76.00	31.87	33.50	550
4	1200	GENDEX GX-CTM (3310) CEILING SUSPENDED X-RAY TUBE CRANE	188.00	108.00	102.00	600
5	320	GENDEX VS-200 SLIM LINE VERTICAL BUCKY STAND	25.00	12.25	84.00	200
6						
7						
8						

GAMMA TECH, INC



129 Wild Berry Lane
Hampstead, NC 28443
Ph (910) 270-9787
Fax (910) 270-9767
www.gammatech.net

© PLANNING DIV.

FILE NAME:
ATC SERIES 3310
REVISION DATE:
07.15.2008