



GENERAL NOTES

- ALL SPECIFICATIONS AND SPACE REQUIREMENTS ARE SHOWN GENERICALLY. CONTACT GAMMA TECH, INC. FOR SPECIFIC REQUIREMENTS PERTAINING TO OPTIONAL EQUIPMENT COMPONENTS AND ADDITIONAL ROOM PREPARATION SPECIFICATIONS.
- MAIN VOLTAGE: 200/208/220/230 VAC NOMINAL (TAP SELECTABLE AT INSTALLATION)
MAINS CURRENT RATING: 35 AMP, WITH DISCONNECT.

EQUIPMENT SCHEDULE

EQUIPMENT IDENTIFICATION NUMBER

	BTU / HOUR	WEIGHT (LBS.)				
		DESCRIPTION	DIMENSIONS (IN.)			
			W	x	D	x
1	2600	LORAD M-IV OPERATOR CONSOLE (WITH SHIELD & OPTIONAL BAR CODE READER)	36.50	17.00	74.50	214
2		LORAD M-IV GANTRY WITH C-ARM	25.50	50.20	84.00	800

GAMMA TECH, INC



129 Wild Berry Lane
Hampstead, NC 28443
Ph (910) 270-9787
Fax (910) 270-9767
www.gammatech.net

© PLANNING DIVISION

FILE NAME:
LORAD MIV

REVISION DATE:
07.17.2008

LORAD

**M-IV
MAMMOGRAPHY SYSTEM**