



# TYPICAL ROOM LAYOUT



SCALE

## GENERAL NOTES

- ALL SPECIFICATIONS AND SPACE REQUIREMENTS ARE SHOWN GENERICALLY. CONTACT GAMMA TECH, INC. FOR SPECIFIC REQUIREMENTS PERTAINING TO OPTIONAL EQUIPMENT COMPONENTS AND ADDITIONAL ROOM PREPARATION SPECIFICATIONS.
- RECOMMENDED CEILING HEIGHT IS 8'-0".
- SUGGESTED AMBIENT ROOM REQUIREMENTS ARE 66° TO 75° F. AT 30% TO 60% R.M. (NON-CONDENSING).
- MAIN INCOMING POWER SUPPLY: 208VAC, SINGLE PHASE, 20AMP, 2-POLE.

## LORAD MULTICARE BREAST BIOPSY SYSTEM SUITE

## EQUIPMENT SCHEDULE

		EQUIPMENT IDENTIFICATION NUMBER				
		BTU / HOUR	WEIGHT (LBS.)			
		DESCRIPTION	DIMENSIONS (IN.)			
			W	D	H	
1	1400	LORAD MULTICARE BREAST BIOPSY TABLE	82.50	66.00	55.00	600
2		LORAD MULTICARE CONTROL CONSOLE	26.50	23.00	59.50	150
3		LORAD DIGITAL SPOT VIEWING AND PROCESSING CONSOLE	34.00	30.00	48.00	50
4	35	LORAD VISIPLEX CCH720 (OR EQUIVALENT) VIDEO IMAGER	17.25	20.75	9.62	38

## GAMMA TECH, INC



129 Wild Berry Lane  
Hampstead, NC 28443  
Ph (910) 270-9787  
Fax (910) 270-9767  
www.gammatech.net

© PLANNING DIV.

FILE NAME:  
LORAD MULTICARE  
REVISION DATE:  
07.17.2008