



TYPICAL ROOM AND EQUIPMENT LAYOUT

GENERAL NOTES

- ALL SPECIFICATIONS AND SPACE REQUIREMENTS ARE SHOWN GENERICALLY. CONTACT GAMMA TECH, INC. FOR SPECIFIC REQUIREMENTS PERTAINING TO OPTIONAL EQUIPMENT COMPONENTS AND ADDITIONAL ROOM PREPARATION SPECIFICATIONS.
- INPUT LINE: 200-240 VAC, (±10%), 35AMP, 50-60HZ
- VERTICAL TRAVEL: 63.5 TO 140 CM (25 TO 55 IN) (MOTORIZED)
 ROTATION: +195° TO -150°
 SID: 65 CM
 PATIENT FACE SHIELD: REMOVABLE

EQUIPMENT SCHEDULE

	EQUIPMENT IDENTIFICATION NUMBER					
	BTU / HOUR	WEIGHT (LBS.)				
		DESCRIPTION	DIMENSIONS (IN.)			
		W	x	D	x	H
1	LORAD SELENIA ACQUISITION WORKSTATION (WITH SHIELD)	38.00	20.50	74.60		
2	LORAD SELENIA GANTRY WITH C-ARM	53.00	50.30	84.00		
3	LORAD SOFTCOPY REVIEW WORKSTATION MONITORS	15.50	21.60	22.60	150	
4	LORAD SOFTCOPY REVIEW WORKSTATION KEYBOARD	18.00	7.00	1.50	2	
5	LORAD SOFTCOPY REVIEW WORKSTATION CPU	8.50	19.25	19.75	55	
6	LORAD SOFTCOPY REVIEW WORKSTATION UPS	6.70	17.30	8.50	53	
7	AGFA DRYSTAR 4500 FILM PRINTER	23.00	29.50	37.00	210	

GAMMA TECH, INC

 129 Wild Berry Lane
 Hampstead, NC 28443
 Ph (910) 270-9787
 Fax (910) 270-9767
 www.gammatech.net

© PLANNING DIVISION
 FILE NAME:
 LORAD SELENIA
 SOFT COPY
 REVISION DATE:
 07.17.2008

LORAD **SELENIA**
 MAMMOGRAPHY SYSTEM