



TYPICAL ROOM AND EQUIPMENT LAYOUT

GENERAL NOTES

- ALL SPECIFICATIONS AND SPACE REQUIREMENTS ARE SHOWN GENERICALLY. CONTACT GAMMA TECH, INC. FOR SPECIFIC REQUIREMENTS PERTAINING TO OPTIONAL EQUIPMENT COMPONENTS AND ADDITIONAL ROOM PREPARATION SPECIFICATIONS.
- INPUT LINE: 200-240 VAC, ($\pm 10\%$), 35AMP, 50-60HZ
- VERTICAL TRAVEL: 63.5 TO 140 CM (25 TO 55 IN)(MOTORIZED)
ROTATION: +195° TO -150°
SID: 65 CM
PATIENT FACE SHIELD: REMOVABLE

EQUIPMENT SCHEDULE

	EQUIPMENT IDENTIFICATION NUMBER	BTU / HOUR		WEIGHT (LBS.)	
		DESCRIPTION		DIMENSIONS (IN.)	
		W	D	x	H
1		LORAD SELENIA ACQUISITION WORKSTAION (WITH SHIELD)	38.00	20.50	74.60
2		LORAD SELENIA GANTRY WITH C-ARM	53.00	50.30	84.00

GAMMA TECH, INC



129 Wild Berry Lane
Hampstead, NC 28443
Ph (910) 270-9787
Fax (910) 270-9767
www.gammatech.net

© PLANNING DIVISION

FILE NAME:
LORAD SELENIA

REVISION DATE:
07.17.2008

LORAD

**SELENIA
MAMMOGRAPHY SYSTEM**