



TYPICAL ROOM LAYOUT



GENERAL NOTES

- ALL SPECIFICATIONS AND SPACE REQUIREMENTS ARE SHOWN GENERICALLY. CONTACT GAMMA TECH, INC. FOR SPECIFIC REQUIREMENTS PERTAINING TO OPTIONAL EQUIPMENT COMPONENTS AND ADDITIONAL ROOM PREPARATION SPECIFICATIONS.
- MINIMUM RECOMMENDED CEILING HEIGHT IS 9'-6".
- SUGGESTED AMBIENT ROOM REQUIREMENTS ARE 66° TO 75° F. AT 30% TO 60% R.M. (NON-CONDENSING).
- MAIN INCOMING POWER SUPPLY: VARIES WITH GENERATOR KW

EQUIPMENT SCHEDULE

EQUIPMENT IDENTIFICATION NUMBER	BTU / HOUR	DESCRIPTION	WEIGHT (LBS.)			
			DIMENSIONS (IN.)			
			W	x	D	x
1	1707	SEDECAL COMPACT SERIES HF GENERATOR CONTROL CONSOLE	17.05	11.73	40.28	75
2		SEDECAL COMPACT SERIES HF GENERATOR LARGE POWER CABINET	23.30	14.17	27.16	165
3	1365	SEDECAL URS_LP LATERAL POSITIONING ARM UNIVERSAL RADIOGRAPHIC SYSTEM	83.07	66.14	104.43	714
4	164	CANON MEDICAL SYSTEMS CXDI-40EG DR SYSTEM 15" FLAT SCREEN MONITOR AND STAND	15.67	5.91	15.55	11
5	250	CANON MEDICAL SYSTEMS CXDI-40EG DR SYSTEM CONTROL PC UNIT (BELOW COUNTER) W/ KEYBOARD & MOUSE	11.81	19.76	21.50	77
6	139	CANON MEDICAL SYSTEMS CXDI-40EG DR SYSTEM - POWER SUPPLY FOR WALLSTAND SENSOR (WALL MOUNTED)	15.35	4.33	5.91	10
7		CANON MEDICAL SYSTEMS CXDI-40EG DR SENSOR UNIT (MOUNTED ON THE WALL)	21.70	21.70	2.70	44

SEDECAL HIGH FREQUENCY URS LP DIGITAL RADIOGRAPHIC SYSTEM SUITE

GAMMA TECH, INC



129 Wild Berry Lane
Hampstead, NC 28443
Ph (910) 270-9787
Fax (910) 270-9767
www.gammatech.net

© PLANNING DIV.

FILE NAME:
URS_LP.DWG

REVISION DATE:
03.12.2010